

NEW

# N-95 Respirator Mask



- Single use particulate respirator mask
- Respiratory protection
- High filtration
- 4-layer construction
- Nose clip on the outside
- Has FDA authorization for emergency use (EUA) during COVID
- Distributed by Mölnlycke and manufactured by Winner Medical, a qualified ISO certified supplier compliant to the Mölnlycke Supplier Code of Ethics

## N-95 Face Mask Test Results Summary<sup>1</sup>

Test Item	Request	Product performance
Carbon Dioxide - Inhalation CO <sub>2</sub> (%)	< 1	0.49-0.59
Breathing Resistance - Inhalation (mbar)	< 0.7 < 2.4	0.32-0.35 1.05-1.12
Breathing Resistance - Exhalation (mbar)	< 3.0	1.80-1.89
Penetration of filter Sodium Chloride Paraffin oil	< 6	0.62-0.79 5.54-5.78
Total inward leakage (%)	≤ 8	Pass

Testing in accordance with BSI COVID-19 filtering face piece technical specification, BS EN 149:2001 + A1:2009

Product No.	Model No.	Description	Case Quantity	Order Quantity (case)	Order Frequency
604-007286	WN-N95FW	Particulate Respirator Mask (N95), Unvalved, folded	20 boxes per case, each box contains 50 masks (1000 masks)		<input type="checkbox"/> Check box to one time order <input type="checkbox"/> 2 months <input type="checkbox"/> 3 months <input type="checkbox"/> Check box to indicate expected no. of orders: (max supply of 6 months from first order date) <input type="checkbox"/> 4 months <input type="checkbox"/> 5 months <input type="checkbox"/> 6 months

Delivery schedule will vary based on production

# About Winner Medical

Winner Medical is a leading manufacturer for disposable wound-care and surgical products in China. Since 1991, Winner Medical has been devoted to manufacturing, R&D and marketing of medical dressings, medical disposables and consumer products. Winner Medical was one of the first companies with full product lines in the disposable medical products market registered with the FDA, as well as numerous ISO certifications. Conforming with pharmacopeia rules in the EU, US, Japan, and China, these products have quickly become the industry benchmark. With multiple local and international honors, awards, and recognitions, Winner Medical has continuously been recognized by the medical device industry. The pursuit of health, sustainable ecology, and higher living standards, enables Winner Medical to continue to capitalize on patented technologies and create great new products.<sup>2</sup>

## IMPORTANT CONSIDERATIONS:

Mölnlycke and Winner Medical are trusted supply chain partners. Winner Medical is in full compliance to the Mölnlycke Code of Ethics, as required by all of Mölnlycke's supply chain partners.

### Quality Standards:

- Meets all of Mölnlycke's internal quality standards for contract manufacturing
- ISO 13485 certified
- Materials:
  - Outer: PP non-woven fabric (50g, white)
  - Inner: Polyester cotton (50 gsm, white)
  - Filter: Fabric ((MA-25X2) 50g, white)
  - Lining: PP non-woven fabric (50 gsm, white)
  - Nose bar: Aluminum
  - Ear loop: Spandex + polyester fiber (white)
  - Green silicone ring
- Color: White
- Country of Origin: China

**References:** 1. Nelson Laboratories test report 1277422-S01.1 2. "About Us" Winner Medical website. Winner Medical. July 21 2020

The Mölnlycke trademarks, names and logo types are registered globally to one or more of the Mölnlycke Health Care Group of Companies. Distributed by Mölnlycke Health Care US, LLC, Norcross, Georgia 30092. © 2020 Mölnlycke Health Care AB. All rights reserved. 1-800-843-8497.

Manufacturer Winner Medical Group Inc.

Lit# 1-0512 7.29.2020

MHC-2020-78844

