

## PFXT SERIES – Class 100 Powder Free Latex Gloves



### PRODUCT OVERVIEW

CTI's PFXT Series are Class 100 Powder Free Gloves. Like all our gloves, they are manufactured for outstanding strength, stretch and fit. They offer a textured surface and are double chlorinated reducing contact sensitizers and protein content. PFXT Series gloves are ISO 5 Clean Room compatible or higher for controlled environments. We offer them in both 10" and 12" long. They are ambidextrous, and have been tested to meet and/or exceed all standards and requirements. Available in sizes Small to X-Large.

### PRODUCT FEATURES

Textured Surface

Low Protein and Contact Sensitizers

Class 100 (ISO 5) Clean Room Compatible or Higher for Controlled Environments

Double Chlorinated

Double Bagged

Tapered, Beaded Cuff

10" and 12" Length, Ambidextrous

Total Protein Content <50ug/g of Glove (ppm) Using Modified Lowry Method

Meets or exceeds FDA requirements with 510(k) registration

Leak Tested

Meets ASTM D3578-05 standards

Sampled in accordance with MIL STD-105E

Complies with U.S. Government 21 CFR 177.2600 for food handling

Qualified for Use in the Pharmaceutical Industry

### SPECIFICATIONS & PACKAGING

Length: **10" and 12" (+/-0.5)**

Thickness: **6.5 mil (+/-1)**

Tensile Strength: **3400 psi/23 mpa**

Elongation: **745%**

Packaging: **Poly Bagged, 100 Gloves/Pack, 1000 Gloves/Case**

Cat No./Sizes:

**10" Sizes:** PFXT101 – Small

PFXT102 – Medium

PFXT103 – Large

PFXT104 – X-Large

**12" Sizes:** PFXT201 – Small

PFXT202 – Medium

PFXT203 – Large

PFXT204 – X-Large