



# Cotton Tipped Applicator

Regular tip, Wood handle



## Description

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. Special construction eliminates chemical migration.

## Dimensions

3.032" Overall Length (77.013 mm)

### Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

### Handle

- 0.100" Dia. (2.54 mm)
- 2.906" Length (73.812 mm)

## Qualities

- Regular cotton tip, no glue adhesive
- Made in USA

## Applications

Critical Environments / Industrial

- Used for cleaning and removing residues

Ref #	Packaging
837-WC NO GLUE Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)




**Puritan®**

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA  
ISO 9001 2008 ISO 13485:2003 CE

Tel 800 321-2313 (US & Canada Only)  
or 207-876-3311

Web [puritanmedproducts.com](http://puritanmedproducts.com)  
Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

### Distributed By:



**HUTCHINS & HUTCHINS™**  
Your Cleanroom Supplier®

Call for more info, pricing & availability | Toll Free: (800) 554-4736 | Local: (540) 949-6663