



HUTCHINS & HUTCHINS
Your Cleanroom Supplier



Products for Every Controlled Environment



TECHNICAL DATA SHEET: **HUTCH-PAK Polypropylene Spunbond Shoe Cover**

Model #	Description	Dimensions (L x W)	Packaged
HHS-28033	Universal	6.93" x 6.5"	3 Easy-Open Perforated Bags 50 Pair / Bag 150 Pair / CS
HHS-28035-2XL	Jumbo (2XL)	18.5" X 7.3"	2 Easy-Pullout Dispenser Boxes 400 Units / CS

DESCRIPTION

Manufactured from Lightweight, Breathable, Non-Woven, Spunbond Polypropylene fabric and 100% Latex-Free. Ideal for transition between 'clean' vs 'unclean' areas. Ultrasonically welded with Anti-Skid Traction and Triple Elastic Runs allow a smooth wrap around the foot helps prevent tripping.

INDUSTRIES / APPLICATIONS

- ◆ Medical / Healthcare Facilities
- ◆ Pharmaceutical Compounding
- ◆ Laboratories
- ◆ Research & Development
- ◆ Food Service
- ◆ Tissue Bank
- ◆ Semi Conductor
- ◆ Biotechnology

BENEFITS & FEATURES

- ◆ Reinforced Corners & Edges
- ◆ Lightweight, Breathable Material
- ◆ Anti-Skid, Seamless Design
- ◆ Wide Opening for Easy Don & Doff
- ◆ Disposable, Single Use
- ◆ High Quality, Low Linting

PHYSICAL CHARACTERISTICS

- ◆ Anti-Static Agent built into Spunbond Polypropylene fabric, Weight: 35 grams/m²
- ◆ Traction Pattern: 100% Latex Free (1.2 grams/m²)
- ◆ Seams: Ultrasonically Welded
- ◆ Manufactured from a single fabric layer, to prevent linting (*no residue*)
- ◆ Color: Light Blue with White Anti-Skid Traction Markings
- ◆ Made in Canada
- ◆ Classification: Class 1 Medical Device

W O M A N - O W N E D S M A L L B U S I N E S S