

#### ■ MAC 10 LEAC Replacement Parts List

ACLE PARTS	PART #
Grommet	63271
Motor 120V	63764-001
Motor 220V	63764-002
Motor 277V	63764-003
Motor Mount Clip	38511-001
Venturi Ring	64115
Blower Wheel	64015
Motor/Blower Assy	
120V	24511-001
220V	24511-002
277V	24511-003
Capacitor	61485
Capacitor Boot	60126
Rocker Switch	63739-002
3 Spd Switch	63775
Foam Pre-filter (20x20)	62981-038
30/30 Pre-filter (20x20)	64125-001

All filter part numbers are based on the standard MAC 10 sizes (2x4 – 23.63x47.63).

If the unit in question is not this size or the part number includes a “Z”, contact the factory for replacement filter information.

#### STANDARD FILTERS PART NUMBERS

- 69600S-001HAPXX FILTER - HEPA 2X4 9” FRAME
- 69600S-001UAPXX FILTER - ULPA 2X4 9” FRAME


#### WELDED RSR GEL SEAL

- 69601-001H FILTER - HEPA 21X45X3.5 GEL SEAL
- 69601-001U FILTER - ULPA 21X45X3.5 GEL SEAL

#### Optional Accessories:

- Fluorescent lighting
- Ionizing bar
- ULPA Filter (Standard and RSR)
- A/C intake collar: 20x20 (10” /254mm)
- 20x20 (12” /304.8mm)
- 20x20 (14” /355.6mm)

#### Distributed By:



**HUTCHINS & HUTCHINS**  
Your Cleanroom Supplier<sup>®</sup>

Call for more info, pricing & availability | Toll Free: (800) 554-4736 | Local: (540) 949-6663