



Validated Sterile SLEEVE/GLOVE SYSTEM 32" LENGTH

TECHNICAL DATA SHEET

PRODUCT INFORMATION

Material	Nitrile Sleeve / Polychloroprene Glove
Color	White Sleeve / Natural Glove
Glove Design	Hand specific hand shape, textured surface
Glove Description	Clean & sterile sleeve/glove system, Nitrile sleeve attached to a size 8.0 or 8.5 hand specific PCP glove by a channel ring and O'ring assembly. Sold by the pair, individually packaged. The out of the bag sterile system is ready to use.
Audit Standards	Manufactured in a facility holding ISO 9001:2015 and PPER 2016/425 Module D
Quality Control	Compliant with glove standards EN420:2003+A1:2009, EN374-1:2016, EN374-5:2016 and EN421:2010. Tested against permeation standard ASTM D6978. Sterilized in accordance with ISO 11137-1. Glove is suitable for use with all foodstuffs. Conforms to Category 3 - Personal Protective Equipment PPE Regulation (EU) 2016/425.
Packaging	One sleeve/glove system per triple bag including size 8 glove (marked L or R) and channel ring/O-ring assembly; 10 systems per white Correx box; Two Correx boxes per outer white Correx box (20 systems)
Storage	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.
Country of Origin	Malaysia
Shelf Life	Three (3) years from date of manufacture.



PHYSICAL PROPERTIES - GLOVE

PROPERTY	TYPICAL VALUES	TESTING METHOD
Freedom from holes	AQL 0.65	EN 374-2:2015
Palm thickness single wall	0.15mm/5.91mil	EN 455-2:2015
Force at break	≥9N	EN 455-2:2015

TESTING INFORMATION - GLOVE

EN ISO 374		
	CHEMICAL	PERFORMANCE LEVEL
K	Sodium Hydroxide 40%	6
O	Ammonium Hydroxide 25%	2
T	Formaldehyde	6


PHYSICAL PROPERTIES - SLEEVE

PROPERTY	TYPICAL VALUES	TESTING METHOD
Cuff thickness single wall	12mil/0.30mm	EN 455-2:2015
Tensile strength	≥14Mpa	ASTM D 6319
Force at break	≥9N	EN 455-2:2015

Tested against ASTM D6978 standard for cytotoxic drugs. Please contact us for results.

ORDERING INFORMATION

Port Size (in/mm)	System Length (in/mm)	Glove SIZE	Product Code
10-12"/254-305mm	32"/813 mm	8	GSG10NIT80
10-12"/254-305mm	32"/813mm	8.5	GSG10NIT85



Distributed By:



Call for more info, pricing & availability | Toll Free: (800) 554-4736 | Local: (540) 949-6663

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates.
© 2019 All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

Warning: No glove provides complete protection against all chemicals. Users must test the gloves against the particular chemicals and environment in which they will be used.

Please see product validation pack or contact Ansell customer service for specific data on use of gloves with cytotoxic drugs. Gloves used for protection against such drugs must be selected specifically for the type of chemicals used.

