

## Laryngoscope Sets

Propper has various sets of laryngoscopes that can assist a medical professional in all situations. Suitable for emergency medical services, professional medical offices or anesthesiology groups; the various set's curved and straight blades will fit most anesthesiologist's needs engaged in adult and pediatric practices. Propper handles are designed to facilitate sterilization and reduce the risk of infection.

**Large Set** - Contains 5 Macintosh blades, 5 Miller blades, 2 Propper blades and 2 handles in a rigid black attaché case.

**Emergency Set** - Contains 1 each of Propper #1, Macintosh #3 and Miller#3 blades, plus a "C" hand "AA" handle in a compact soft case.

**Flip-Tip® Set** - 2 Flip-Tip® blades and 1 handle in a rigid compact attaché case.

All sets are available with either Superlume® LED or Xenon bulbs.

Reorder no.	Description
199176	Large
199149	Large LED
199168	Flip-Tip®
199171	Flip-Tip®
199172	Flip-Tip® LED
199173	Flip-Tip® LED
199175	Emergency
199178	Emergency LED



## Meeting The Goal

Propper's goal throughout its prestigious history has always been to provide the medical professional with safe, reliable and effective instruments. No words can convey our commitments more strongly than holding and using our laryngoscope products. The superior quality and flawless craftsmanship will be immediately evident. The large number of medical facilities throughout the world currently using Propper laryngoscope products confirms that our commitment to the goal has been met time and time again.



ARE YOU USING THE PROPPER LARYNGOSCOPES?

Offering an extensive family of laryngoscope handles, blades & sets for Anesthesiologists, Emergency Medical Technicians & other medical specialists.

## Laryngoscope Handles

Propper Integrated Fiber Optic Laryngoscopes have been considered the quality standard in the operating room since their introduction. The 3 handle sizes, "C", "Stubby" & "AA" are available with xenon bulbs or the new Superlume® LED light sources. Propper laryngoscopes provide you with one brilliant advantage after another: brighter illumination, ease of cleaning & sterilization and longer lifetime. "C" & "AA" handles feature a long life patented switch mechanism and a removable inner compartment allowing hassle-free cleaning and sterilization. All 3 handles fit any fiber optic blade compliant with the "green system" standard.

The state of the art Superlume® LED handles have the brightest output of any current handle. Together with its high color temperature, the effect creates the best conditions for intubation. Moreover, LED lifetimes will make bulb replacement a thing of the past.

US Patent #5,542,905

Reorder No.	Description
199581	Stubby LED
199578	Stubby Xenon
199582	C LED
199150	C Xenon
199583	AA LED
199152	AA Xenon



Stubby



AA handle



C handle



## Laryngoscope Blades



Mac



Miller



Propper



Flip-Tip®

- Fully compliant with infection control requirements
- Integrated fiber optic streamlined design
- High performance fiber optic bundle
- Large diameter fiberoptics with unsurpassed light transmission
- Superior material and surface treatment
- Designed to be sterilized using steam, Steris®/Sterrad®, or EO
- Extended lifetime insures lowest cost per use
- Virtually maintenance free
- Reinforced tip
- Flawless pinhole free Ni/Cr plating
- Select quality 304 stainless steel
- Epoxy sealed channel heat resistant up to 392°F
- Patented integrated design
- Reinforced ball lock system
- Integrated design minimizes bodily fluid entrapment.
- Fits any "Green System" handle.

Reorder No.	Description	Size
199153	Propper Blade	0
199147	Propper Blade	1
199161	Macintosh Blade	1
199144	Macintosh Blade	2
199145	Macintosh Blade	3
199146	Macintosh Blade	4
199162	Macintosh Blade	5
199148	Miller Blade	00
199156	Miller Blade	0
199157	Miller Blade	1
199155	Miller Blade	1.5
199158	Miller Blade	2
199159	Miller Blade	3
199160	Miller Blade	4
199170	Flip-Tip Blade	2
199166	Flip-Tip Blade	3
199167	Flip-Tip Blade	4

## Flip-Tip® For Difficult Airways

The Propper Flip-Tip® blade is a modification of the traditional Macintosh blade designed to improve visualization of the larynx, and can assist in both non-complicated and difficult intubations. The instrument is especially suitable in situations of unpredictable airway conditions in emergency rooms, ambulance service, etc.

Manipulation with the Propper Flip-Tip® laryngoscope does not require additional training and reduces the risk of teeth breaking and laryngeal tissue injury. The Flip-Tip® blade can simplify elevation of the epiglottis in cases of spine or head trauma, immobilized neck, reduced thyromental distance, swollen tissue of the larynx, obesity, enlarged tongue, prominent upper teeth and poor mouth opening.



visit the new [www.proppermfg.com](http://www.proppermfg.com)

email [customerservice@proppermfg.com](mailto:customerservice@proppermfg.com)